

EMR Connect for the Tablo® Hemodialysis System

Say “Goodbye” to manual charting



Nurses spend up to

41%

of their time on documentation.*

Streamline your dialysis workflows from charting to billing.

- ✓ Automated patient charting and documentation for Tablo to save staff time
- ✓ Real-time data transfer from Tablo to your EMR platform to eliminate manual charting
- ✓ Accurate flowsheets in patient chart designed to minimize documentation errors

*De Groot K, Munster AM, et al. Nursing documentation and its relationship with perceived nursing workload: a mixed-methods study among community nurses, MBC Nursing, (2022) 21:34.

Streamlined connection from the patient to your EMR.



Patient treatment



Secure, cloud-based connection



Real-time patient charting in EMR

Technical overview

DATA OPTIONS

- 70+ treatment fields for charting
- Data transmission timing defined by customer
- Event-based data and alarms sent in real-time
- One-way dialysis treatment data transmission

CONNECTIVITY

- HL7v2 and API/JSON communication options
- Secure, dedicated, cloud-based connection
- Support for HIPAA compliance

SUPPORTED EMR VENDORS

- Major vendors such as Epic, Meditech, Gaia and Cerner
- Full vendor list upon request



Want to learn more about
EMR Connect for Tablo?

Contact us at 844-MYTABLO,
cs@outsetmedical.com or
www.outsetmedical.com/contact-us